The objective to reach with this uncomplicated module is the study in depth of the **equation of Bernoulli** and its demonstration.

This apparatus is based on Bernoulli's principle, which describes the behaviour of a laminar flow moving throughout a conduit and considers that in an ideal situation, for a fluid (without viscosity nor friction) in a regime of circulation through a closed conduit, the energy remains constant throughout its route.

The equipment comes with a Venturi tube, implemented to carry out experiments for the demonstration of the theory based on the theorem of the conservation of mechanical energy.

The regulation valve allows to work with a variety of flow rates giving rise to different scales in differential pressure.

The pressure results are shown in the multi-tube manometers installed in the equipment. Therefore, the readings of pressure along the conduit can be obtained quite easily.


The measure of flow rates is carried out thanks to the volumetric tank in the hydraulic bench (required), which also enables the study of the relation between the pressure drop and the speed of the fluid.


Venturi is easy to install and fully transparent.

Fasy connection to hydraulics bench.

8 tapping Multimanometer.

Pitot to know the speed of the fluid at each point.

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9

The manual shows clearly and with a lot of images, the hole process to operate the equipment.

**LEARNING OBJECTIVES**

Some of the learning objectives of this equipment are the following:

- Demonstration of the equation of Bernoulli in a Venturi nozzle.
- Calculation of the losses in a Venturi nozzle.
- Study of the static, dynamic and total pressures.
- Study of the Venturi nozzle as a flow measurer and calculation of the discharge coefficient.

**TECHNICAL DATA****Inner diameters:**

- Main pipe  $\varnothing = 25$  mm.

**Venturi nozzle:**

- Min. diameter  $\varnothing 10$ mm.
- Inlet angle  $14^\circ$ .
- Outlet angle  $21^\circ$ .

**Pressure measurement points:**

- a =25 mm
- b =20 mm
- c =10 mm
- d =14 mm
- e =18 mm
- f =20 mm
- g =25 mm

**Pressure gauges:**

- 8 tube manometer, measurement range 300 mmWC.

**REQUIREMENTS**

- DIKOIN hydraulic bench.