With the HD 10.1 equipment, it is a question of knowing the water retention capacity of a soil using the hydrological balance, differentiating each one of its elements.

By means of the simulation of a rain on a soil, the following elements of the hydrological cycle are distinguished:

- Surface runoff: Once the soil is flooded, surface runoff begins, collecting it externally and proceeding to its measurement by evaluating its quantity in volume.
- Subsurface and subterranean runoff: The infiltrated water is collected by the bottom of the container, measuring this runoff by its volume.
- The rest of the water is absorbed by the soil that can be known by performing a weighing of the same before the rain and another after. The difference corresponds to the volume of stored water.

By establishing the equality of volumes, the values of surface runoff, ground runoff and soil runoff are known.

The equipment has a digital balance on which the floor model is set, so that the weight variation can be continuously displayed.

The equipment allows the height of the diffuser to be adjusted to position it correctly according to the flow rate used.

Two drums are supplied to collect and measure the amount of surface and groundwater runoff.

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The manual shows clearly and with a lot of images, the hole process to operate the equipment.

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The instruction manual explains and shows all the theoretical foundations, as well as all the mathematic expressions used during the experimentation.

17

DIKOIN Ingeniería S.L.U.

+034 946 55 15 35 / www.dikoin.com / info@dikoin.com
Dikoin reserves the right to make changes in technical and product specifications without prior notice.

LEARNING OBJECTIVES

- Separation of the components of the hydrological cycle by volume or height.
- Soil retention capacity.
- For a short duration of experiments, evaporations are neglected but this duration can be prolonged and obtatin the soil evaporation (a thermometer is required to carry out this experiment).
- For a long duration of experiment placing a plant or crop, you can determine the evapotranspiration and the point of wilt.

TECHNICAL DATA

• Peripheral drive pump (which simulates precipitated rain volume on the floor model):

Maximum flow: 10 l/minH max: 42 m a 50HzPower absorbed 0,25CV

- Tanks
 - Test container 501
 - Stackable tank 12l
 - Drum 20I
- Flowmeter (to measure the discharged water): 16-160l/h
- Surface excess flow rate measurement system.
- Infiltrated water measurement system.
- Digital scale 60Kgx2g

REQUIREMENTS

• Input: 230V/50Hz.