The TC 01.1 equipment is the core of the whole heat exchanger. TC 01.1 is the module that provides hot and cold water to the heat exchangers, in addition to measuring the temperatures and flow rates for each element.

All device connections are fast, self-sealant connectors that allow for a quick and simple change of exchangers without any loss of fluid. The connections for hot and cold water are clearly differentiated to prevent mistakes.

The equipment unit has a tank for hot water with 4,5 litre capacity, as well as electronic controllers both of temperature and water level. An electromechanical valve fills up the tank automatically as needed. The water storage system is protected against overheating, low water level and overflowing.

The pumping system has a bypass line, which makes the pump operation smoother and the conditions in the tank steadier.

TC 01.1 HEAT EXCHANGER SUPPLY

HOT FLOW

T supply 13.9 °C

T return 13.5 °C

Flow 0.0 l/min

Cold Flow

T supply 14.1 °C

T return 19.7 °C

Flow 0.0 l/min

SHELL & TUBE HEAT EXCHANGER

JACKETED VESSEL WITH STIRRER & COIL

TUBULAR HEAT EXCHANGER

PLATE HEAT EXCHANGER

CALIBRATE

EXIT

The equipment includes a PC with software to manage the equipment. In it, all the control points of the equipment are shown, and the data acquisition is permitted in manual or automatic mode.

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The equipment has an automatic sensor calibration system.

Optional Accessory: TC 01.2 - PLATE HEAT EXCHANGER

In the plate heat exchanger, the hot and cold flows alternating sides pass through the gaps left by the plates, thus resulting in heat transfer.

The advantage of this type of heat exchanger is a compact configuration, and therefore are suitably used in confined spaces.

The plates have a geometry that causes a turbulence in the fluid, improving heat transfer.

Optional Accessory: TC 01.3 - SHELL TUBE HEAT EXCHANGER

The beam exchanger tubes is one of the most used in industry exchangers.

In this exchanger, the cold fluid passes through a series of parallel tubes grouped, and the heated fluid through the chamber containing small, thus resulting in heat transfer.

The advantage of this type of heat exchanger is a compact design and the ability to work at higher pressures than other designs.

This exchanger can operate with co-current or countercurrent flows.


Optional Accessory: TC 01.4 - TUBULAR HEAT EXCHANGER
The heat exchanger of concentric tubes is simpler in design.

There are two parallel tubes through which cold fluid passes, inside which there is another pipe of smaller diameter by passing the heated fluid, thereby producing heat transfer.
The advantage of this exchanger is its simple design.

The exchanger is arranged in two halves, and has incorporated thermocouples at midpoints, so as to significantly improve the learning experiment, because you can clearly see the change in temperature over the heat exchanger.

This exchanger can operate with co-current or countercurrent flows.

Optional Accessory: TC 01.5 - DOUBLE JACKETED VESSEL AND COIL HEAT EXCHANGER

This type of exchanger is usually used in the quimic and process industry, when a very well adjusted temperature is needed.

The exchanger can work with the vessel, or with the coil, and also there is the possibility to work with a continuous flow in the vessel, os heat a defined quantity of water.

The exchanger has also a temperature sensor which measures continuously the fluid temperature into the vessel, as well as a variable speed mixer, to study the differences on heat exchange.

TC 01.1 - HEAT EXCHANGER SUPPLY UNIT

LEARNING OBJECTIVES

- Demonstration of heat transfer.
- Comparison of different types of heat exchangers.
- Comparison of results with flows co-current and countercurrent.
- Transfer coefficient measurement, the effects of flow rate and temperature differential.
- Calculation of energy balance and efficiency.

TECHNICAL DATA

- Adjustable from 0 to 1.5 kW heater from the computer.
- Peripheral Pump:
 - Maximum flow: 10 l / min (5m.c.a.)
 - Power consumption: 180W
- Hot water circuit with Bypass.
- Hot water maximum temperature: 60°C.
- Maximum hot water flow rate: 5 l / min
- The unit is supplied with an electronic and computerized control system and representation, including computer.
- Dimensions: 730x1100x530 mm

ACCESSORIES

- TC 01.2 PLATE HEAT EXCHANGER.
- TC 01.3 HEAT EXCHANGER TUBE BEAM.
- TC 01.4 HEAT EXCHANGER TUBES IN CONCENTRIC.
- TC 01.5 DOUBLE COIL HEAT SHIRT AND COIL.

Note:

The heat exchangers are not included in the 01.1 TC team.
01.1 The TC team needs at least one exchanger to operate.
(The plate heat exchanger shown in the image is not included with the unit)

REQUIREMENTS

- Electrical connection 230V/50Hz
- Water input minimum of: 5 l/min
- Waste water connection